



CATACHEM

FRUCTOSAMINE

PRODUCT / SERVICE INFORMATION

Catachem, Inc. Introduces the Fructosamine Single Working reagent(No Preparation Required) In-Vitro Diagnostic (IVD) Chemistry Reagent Kit.

Catachem Fructosamine two-reagent system has been enhanced to make these reagents available in a single ready-to-use working reagent without the need for preparation. This new Fructosamine reagent packaging eliminates transferring of Reagent-2 into Reagent-1 (old system). The new reagent has a shelf stability of over two years and working stability of six months, upon opening of container and stored at 2-8° C.

The new Fructosamine reagent packaging method is in direct agreement and works with the same instrument parameters currently in use for Catachem Fructosamine two-reagent system.

For technical assistance please call Catachem, Inc. 203-335-9277 or Fax 203-335-9977.